



Ready-to-use Surface Cleaner & Intermediate Level Disinfectant, General Virucide, Bactericide, Tuberculocide.

One-step cleaner disinfectant formulated with proprietary hydrogen peroxide technology to deliver faster, more effective cleaning performance.

Meets Health Canada's standard for cleaning and decontaminating objects and surfaces soiled with blood and body fluids.

Broad Spectrum Disinfection Action

BACTERICIDAL and VIRUCIDAL in one minute.
TUBERCULOCIDAL in one minute.

Improved Cleaning

The patented technology** contains a high performance surfactant system that delivers faster and better cleaning results. No rinsing is required in most applications for improved productivity.

Environmentally Responsible

No VOC's (Volatile Organic Compounds), no added fragrance, no NPE's (nonylphenol ethoxylates). The active ingredient breaks down to water and oxygen.

Non-Corrosive Formula

Can be used to clean most hard, non-porous environmental surfaces with confidence. User Safety—non-irritating to skin and eyes and contains no substances which at the given concentration, are considered to be hazardous to health per Health Canada standards.



Johnson wax
PROFESSIONAL



OXIVIR™ Tb

General Virucide, Bactericide, Tuberculocide, Sanitizer



Profile

WHAT IS IT?	A revolutionary, proprietary one-step cleaner disinfectant for spray and wipe cleaning applications.
WHAT DO I USE IT ON?	Use daily on the following hard, non-porous environmental surfaces: vinyl, plastic surfaces, glazed ceramic, glazed porcelain, chrome, laminated surfaces (and baked enamel surfaces associated with floors), walls, ceilings, tables, chairs, countertops, telephones, fixtures, glazed tile, toilet exteriors, urinal exteriors, sinks found in health care facilities, hospitals, patient rooms, operating rooms, shower rooms, and locker rooms. It can be used to pre-clean and disinfect hospital items: wheelchairs, hospital patient bedrails and linings, patient bed linings, vinyl mattress covers, wash basins, bed pans, non-critical medical equipment surfaces – almost any washable food and non-food contact surface where disinfection is required. A potable water rinse is required when disinfecting food contact surfaces. Do not use on glasses, dishes and utensils.
How do I use it?	Pre-clean heavily soiled areas. Apply solution by cloth or disposable wipe to hard, non-porous environmental surfaces. All surfaces must remain wet for 1 minute. Wipe surfaces dry, or rinse.

Specifications

Active Ingredient: Hydrogen Peroxide 0.5%; **DIN Reg. No.:** 02283530 (Wipes), 02283522 (Liquid)
Color / Form: Clear colorless liquid; **pH:** 2.5 – 3.5; **Fragrance:** No fragrance
Cleaning Performance: Excellent; **Phosphates:** None; **Specific Gravity:** (water = 1) 1.01

A POWERFUL DISINFECTANT AGAINST:

VIRUCIDAL: Proven effectiveness against the Poliovirus Type 1, Strain Brunhilde (ATCC VR-1000) which allows for a general virucidal claim against most enveloped and non-enveloped viruses, HIV-1 Human Immunodeficiency virus (HIV), Strain HTLV-III B (HIV-1), Feline Calicivirus (ATCC VR-782), as a surrogate for Norwalk and Norwalk-like viruses • Human Coronavirus (ATCC VR-740) • Influenza A, Strain Hong Kong (ATCC VR-544) • Herpes Simplex Virus, Type 1 (HSV-1) (ATCC VR-733) • Herpes Simplex Virus, Type 2 (HSV-2) (ATCC VR-734) • Rhinovirus Type 37, Strain 151-1 (ATCC VR-1147) • Rotavirus, Strain WA • Hepatitis C Virus (HCV)

BACTERICIDAL: *Pseudomonas aeruginosa* (ATCC 15442)

• *Staphylococcus aureus* (ATCC 6538) • *Salmonella choleraesuis* (ATCC 10708)
 • *Escherichia coli* 0157:H7 (ATCC 35150) • *Staphylococcus aureus* MRSA (ATCC 33592)
 • *Gentamycin and Methicillin resistant (MRSA)* • *Enterococcus faecalis* VRE (ATCC 51575)
 • *Acinetobacter baumannii* (ATCC 19606)

BROAD-SPECTRUM SANITIZING: *Klebsiella pneumoniae* (ATCC 4352) • *Salmonella choleraesuis* (ATCC 10708) • *Pseudomonas aeruginosa* (ATCC 15442) • *Staphylococcus aureus* (ATCC 6538)
 • *Staphylococcus aureus* MRSA (ATCC 33592) • *Enterococcus faecalis* VRE (ATCC 51575)
 • *Escherichia coli* 0157:H7 (ATCC 35150)

TUBERCULOCIDAL: *Mycobacterium terrae* (ATCC 15755)

Safety Reminder

Please make sure your employees read and understand the product label and Material Safety Data Sheet before using this product. The label contains directions for use, and both the label and MSDS contain hazard warnings, precautionary statements and first aid procedures. MSDS's are available online at www.johnsondiversey.com or by calling 1-888-352-2249.

HMIS/NFPA HAZARD RATING

Health 0
 Flammability 0
 Reactivity 0

Extra Large Wipes!
 20 cm x 23 cm



Packaging

PACK / SIZE	PRODUCT CODE
RTU - 12 / 946 mL RTU Bottles (includes two trigger sprayers per case)	4277293
Wipes - 12/160 ct. Canisters Wipe Dimensions: 20 cm x 23 cm	4444151

Johnson wax
 PROFESSIONAL

For further information on special applications or other Johnson Wax Professional products, contact your JohnsonDiversey representative at 1-800-558-2332 or visit our web-site at www.johnsondiversey.com.